

Thin Film Measurement Software TF ProVis Pro V3.0

The *TF ProVis Pro* is an easy-to-use software application for fast and precise in-line thin film analysis under Windows XP and Vista. Single or 3-axes traverse systems are supported as well as multiple digital and analog ports.

The software uses an improved Fast-Fourier Transformation (FFT) algorithm to determine the film thickness from white-light interference spectra created by thin transparent layers. Layer thickness from below $1\mu\text{m}$ up to nearly $200\mu\text{m}$ can be determined in conjunction with tec5 *MultiSpec* spectrometer systems. The film thickness is computed in real-time and is displayed in various on-line charts. The results are logged to an ASCII file during the measurement. Single as well as double layer structures can be examined.




Fig. 1: Screenshot TFProVis Pro

This software has been created by a partner company with long-time expertise in thin film analysis. It is compatible to tec5 MultiSpec desktop USB spectrometer systems as well as OEM USB operating electronics.

Features and Specifications

Film Thickness Determination Features

- Range: approx. 1 to 150 μ m optical thickness - depends on used spectrometer module respectively the useful spectral range and the optical resolution
- Evaluation of interference spectra by Fast-Fourier Transformation (FFT)
- Run-time-optimized algorithm, evaluation time less than 1ms
- Special algorithm for highly accurate sub-pixel determination of FFT peak position
- Selectable spectral range for interference spectrum evaluation
- Consideration of refractive index and its dispersion (Cauchy dispersion function)
- Selectable film thickness evaluation range within the FFT spectrum for fully automatic determination of double layers

Options for Measurement and Visualization

- Support of triple axes scanning bridges to perform x/y/z coordinate controlled, fully automated measurements of rectangular thickness profiles with up to 250 x 250 values per scan
- Support of single axis traverses for continuously scanning forwards/backwards in cross-direction of a web coater
- Support of 8-channel TTL or 4-channel analog outputs for reporting film thickness values and high/low limits
- Manually or fully automatic triggered spectral data acquisition and film thickness determination, trigger by timer, TTL or traverse controlled
- Real-time representation of interference and FFT spectrum during measurement
- Real-time representation of the film thickness results as trend graphics and bar charts
- Film thickness results are logged to read-shared text files, accessible by third-party software during measurement
- Logs up to 62,500 thickness results to read-shared text files
- Logs up to 62,500 spectra as 'Spectra-Recorder', permits a subsequent off-line film thickness re-evaluation with changed evaluation parameters
- Saves all parameter settings to individual parameter files
- Password protection of parameter setups and special user access rights for each document individually
- Quick access to last used parameters and 'Spectra-Recorder' files

Software Features

- Multi-threaded and Multiple Document Interface handling
- Shell registration for drag-and-drop operation of the document files
- Minimal requirements for resources and memory
- Programmed in Visual C++ by use of the Microsoft Foundation Classes (MFC)
- Consideration of the Microsoft Application Design Guide: menu toolbar, status bar, tool tips, HTML on-line help
- fully supported Windows XP Themes and multi-monitor use
- Software documentation as detailed user User's Manual and Quick Reference Guide

Minimum Hardware and Software Requirements

- PC with at least Pentium 4 processor
- Windows XP or Vista operating system (NT not supported)
- 128 MB RAM, 256 MB or more recommended
- CD ROM drive for installation
- Graphics adapter with 1024 x 768 pixels (1280 x 1024 recommended)
- tec5 MultiSpec Pro spectrometer system with USB high-speed interface