

New Releases tec5 Software Tools

There are new releases for the following software products:

Order ID	Product Designation	Version
11-2010022-00	SDACQ32MP SDK	3.0.1
11-2010029-00	SDPROC32 SDK	2.0.1
11-2010023-00	SDACQ32MP LabVIEW SDK	3.0.1
11-2010026-00	AdminTool	2.0.2
11-2010100-00	MultiSpec Pro	3.5.02

Main new features are the implementation of the Interface Electronics PD-ETH01 and COE-USB11. Some Interface Electronics and electronic components quite old are no longer supported:

- PD-ISA16V3, PD-10401V1, MOE-V1, DPU-USB / MMS-CCD
- FEE-003, FEE-HR, FEE-NIR, FEE-CCD
- MUX-4A

Windows XP / Vista / Windows 7 are supported. The 32 Bit as well as the 64 Bit version of Windows 7 are supported. However, the software tools are realized as 32-Bit software.

MultiSpec Pro software does not run under Win7 yet, this feature will follow with version 3.6. (release planned for end of Q2, 2011).

We recommend to update your software tools.

Please visit <http://www.tec5.com/engl/support/software/index.html> for downloads or contact tec5USA for sending a new InfoCD.